

Dermatology

Derma Roller



We offer easy to operate and maintain Skin Derma Roller or Micro needle roller system to the patrons, which is widely appreciated these days. The micro needle skin nursing system is manufactured using high quality raw materials procured from the most reliable vendors of the industry. They are quality tested several times before delivery to the customers and are available to our customers at very reasonable rates.

General Specifications

- Material of needle: Titanium / SS depending on the model
- Material of handle & roller: PC+ABS
- Quantity of needle: 200 Nos. / 540 Nos depending on the model
- High sealing sterilization packing

Applications of derma roller

- Acne scar 0.5 mm, 0.75 mm, 1.0 mm, 1.5 mm
- Anti age & anti wrinkle 0.5 mm, 0.75 mm, 1.0 mm, 1.5 mm
- Stretch marks 1.0 mm, 1.5 mm
- Hyper pigmentation 1.5 mm, 2.0 mm
- Acne scar deep 1.0 mm, 1.5 mm, 2.0 mm, 2.5 mm
- Hair loss & Regrowth 0.5 mm
- Cellulite 1.0 mm, 1.5 mm, 2.0 mm, 2.5 mm
- Cell regeneration 0.25 mm, 0.5 mm, 0.75 mm, 1.0 mm, 1.5 mm, 2.0 mm, 2.5 mm

Advantages of micro-needling

- Short healing period
- Performed under topical anaesthesia, if required
- No permanent skin damage
- Thin or pre lazed skin can be treated
- Minimal risks, safe and no risk of burns
- Suitable for all type of skin
- All areas of the face, neck and body can be treated
- Cost is less than laser treatments

