

Cosmetology

Hair Transplantation Chair

Hair Transplantation chair is an AC-powered and intended to position a patient to perform routine procedures in most of the specialties of healthcare. The chair is motorized. This chair uses a low voltage motor-driven linear actuator and actuated by a handheld controller to change the position of the seat from sitting to sleeping position. Procedure Chair increases the patient comfort considerably and therefore help for successful procedures.

Available with

- Two motor technology to control up + Down movements and synchronous backrest + Footrest movements using wired remote controller.
- Three motor technology to control up + Down movements, backrest & Footrest movements using wired remote controller
- Medical grade foam and fabric
- Armrests for patients in sitting position

HT 200

Optional Detachable headrest, Optional forehead rest with tilting mechanism, Hand rest for prone position

Optional Features for HT 300:

- ◆ Halogen Procedure Light
- ◆ VGA Camera
- ◆ DVR for storing the procedure
- ◆ 7 inch LED Display
- ◆ Instrument tray

